



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mounting Frame

**Product Code :**  
LBNY-0005-101000116



**Description :**

Frame for mounting of experiments in statics, strength of materials and dynamics.  
Stable on laboratory desktops or workbenches.  
Frame supplied disassembled.  
Sturdy sectional steel double frame, welded.  
Easy, exact mounting of all components by precision clamp fixings.

**Technical Specification :**

Mounting Frame



**Laboratory Instrument India**