



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacuum Pump

Product Code :
LBNY-0002-690006



Description :

Vacuum Pump

Technical Specification :

Integrated pump body structure, 15 micron of ultimate vacuum.
Innovative appearance design with increased area for heat emanation to realize low temperature increase while in running.

Big sight glass, low oil level and forced-feed oil design to avoid running without oil and ensure reliable use..

Vacuum Pump, Vacuum Pump Equipments, Vacuum Pump Tools, Vacuum Pump Tool Kits, Vacuum Pump Manufacturers, Vacuum Pump Suppliers from India, China, Kenya



Laboratory Instrument India