



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Parabolic Arch

**Product Code :**  
LBNY-0005-101000118



**Description :**

Dial gauges record the deformation of the arch under load  
Four sets of weights to compensate for the reactions of a fixed support  
Loading of the arch with a distributed load by way of 7 evenly distributed weights or by point loads  
Investigation of a parabolic arch, optionally statically determinate 1 fixed support, 1 movable support or indeterminate two fixed supports.  
Experimental set-up in frame.  
Storage system to house the components.

**Technical Specification :**

Parabolic Arch



**Laboratory Instrument India**