



Email : sales@laboratoryinstrumentindia.com

Product Name :
Parabolic Arch

Product Code :
LBNY-0005-101000118



Description :

Dial gauges record the deformation of the arch under load
Four sets of weights to compensate for the reactions of a fixed support
Loading of the arch with a distributed load by way of 7 evenly distributed weights or by point loads
Investigation of a parabolic arch, optionally statically determinate 1 fixed support, 1 movable support or indeterminate two fixed supports.
Experimental set-up in frame.
Storage system to house the components.

Technical Specification :

Parabolic Arch



Laboratory Instrument India