



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shear Force Diagram

Product Code :
LBNY-0005-101000121



Description :

Measurement of shear force in beam by low-friction hinge with one degree of freedom.
Investigation of shear force on beam mounted on two supports.
Position of hinge at $1/3$ span.
Loading of beam by one to three point loads.
Two bearing supports.
Force gauge to indicate shear force.
Adjuster nut for horizontal alignment of beam.
Storage system to house the components.

Technical Specification :

Shear Force Diagram



Laboratory Instrument India