



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Gerber Beam

**Product Code :**  
LBNY-0005-101000122



**Description :**

Beam subjected to moving load or loaded down by single weights.  
Gerber beam comprising two cantilever beams and one suspended beam.  
Suspended beam mounted by articulated supports on cantilever beams.  
Investigation of the lines of influence in a statically determinate Gerber beam under different loading.  
Four movable supports with force gauge.  
Experimental set-up in frame.  
Storage system to house the components.

**Technical Specification :**

Gerber Beam



**Laboratory Instrument India**