




Email : sales@laboratoryinstrumentindia.com

Product Name :
Cable Under Dead Weight

Product Code :
LBNY-0005-101000123



Description :

Height of a chain wheel adjustable for unsymmetrical experimental set up.
Cross-arm with scale to hold chain wheels and rule to measure vertical sag of chain.
Graduated weight sets.
Determination of the catenary of a free hanging cable.
Adjustable chain wheel axle spacing.
Storage system to house the components.
Experimental set up in frame.
Symmetrical and unsymmetrical experimental set-up possible.
Roller chain as cable with 2 ball bearing-mounted chain wheels.

Technical Specification :

Cable Under Dead Weight



Laboratory Instrument India