




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cable Under Dead Weight

**Product Code :**  
LBNY-0005-101000123



**Description :**

Height of a chain wheel adjustable for unsymmetrical experimental set up.  
Cross-arm with scale to hold chain wheels and rule to measure vertical sag of chain.  
Graduated weight sets.  
Determination of the catenary of a free hanging cable.  
Adjustable chain wheel axle spacing.  
Storage system to house the components.  
Experimental set up in frame.  
Symmetrical and unsymmetrical experimental set-up possible.  
Roller chain as cable with 2 ball bearing-mounted chain wheels.

**Technical Specification :**

Cable Under Dead Weight



**Laboratory Instrument India**