



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Cooler (Cold and Hot)

Product Code :
LBNY-0002-690007



Description :

Water Cooler (Cold and Hot)

Technical Specification :

This product is installed at homes, offices, hotels and restaurants for providing safe drinking water. The dispenser is capable to offer both hot & cold water.

Water Cooler (Cold and Hot), Water Cooler (Cold and Hot) Equipments, Water Cooler (Cold and Hot) Tools, Water Cooler (Cold and Hot) Tool Kits, Water Cooler (Cold and Hot) Manufacturers, Water Cooler (Cold and Hot) Suppliers from India, China, Kenya



Laboratory Instrument India