



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Shear Force and Bending Moment Diagrams

**Product Code :**  
LBNY-0005-101000129



**Description :**

Measurement of shear force and bending moment in beam by low-friction hinge with two degrees of freedom  
Determination of shear force and bending moment on beam mounted on two supports  
Force gauges to indicate shear force and determine bending moment  
Bending moment determined by force measurement and lever arm  
Position of hinge at 1/3 span  
Two bearing supports  
Loading of beam by one to three point loads  
Steel rule to determine positions of point loads  
Adjuster nuts for horizontal alignment of beam  
Storage system to house the components.

**Technical Specification :**

Shear Force and Bending Moment Diagrams



**Laboratory Instrument India**

