



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shear Force and Bending Moment Diagrams

Product Code :
LBNY-0005-101000129



Description :

Measurement of shear force and bending moment in beam by low-friction hinge with two degrees of freedom
Determination of shear force and bending moment on beam mounted on two supports
Force gauges to indicate shear force and determine bending moment
Bending moment determined by force measurement and lever arm
Position of hinge at 1/3 span
Two bearing supports
Loading of beam by one to three point loads
Steel rule to determine positions of point loads
Adjuster nuts for horizontal alignment of beam
Storage system to house the components.

Technical Specification :

Shear Force and Bending Moment Diagrams



Laboratory Instrument India

