



Email : sales@laboratoryinstrumentindia.com

Product Name :
Forces in a Crane Jib

Product Code :
LBNY-0005-101000132



Description :

Tensile and compressive forces in a planar central force system based on the example of a crane jib.
Various jib forms possible.
Loading with weights set, up to 50N.
Steel weights, surfaces galvanised.
Handles to aid transportation.
Stainless steel retaining bar.
Integrated spring balances in the bars.
Max. load on crane jib 50N.
Sturdy metal frame.
Box to house the components.

Technical Specification :

Forces in a Crane Jib



Laboratory Instrument India

