




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Water To Air Heat Transfer Module

**Product Code :**  
LBNY-0005-10200031



**Description :**

A Cooling fan is supplied with a speed controller to control the airflow over the heat exchanger. Water temperature is also measured via the inbuilt thermocouples of the base unit. Thermocouples are placed before and after the fan to measure air temperature.

**Technical Specification :**

Water To Air Heat Transfer Module



**Laboratory Instrument India**