



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concentric Tube Heat Exchanger Apparatus

Product Code :
LBNY-0005-10200036



Description :

Self sealing quick release connections enable rapid connection to and conversion from parallel to counter current flow. An extended version with 3 pairs of intermediate points and a total of 10 thermocouples is also available.

A clear acrylic tube containing cold water surrounds an inner stainless steel tube in which the hot water flows. Six thermocouples measure hot and cold inlet, mid-point and exit temperatures. Total heat transfer area of approximately 24000mm².

Technical Specification :

Concentric Tube Heat Exchanger Apparatus



Laboratory Instrument India