



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydrodynamic Channel

Product Code :
LBNY-0005-1030004



Description :

Included accessories:

Pitot tube and level meter.
Thin edged wall weir without contraction.
Four-column pressure gauge.
Radial gate.
Vertical gate.

Technical Specification :

Technical Specification:

Gauges:
Multi-column pressure gauge for pressure readings located beneath the channel.
Flow meter:
Electronic flow meter with display.
Working section
Length: 2500mm
Height/ width: 300/80mm
Centrifugal water pump
Max. flow: 24m³/h
Power: 1.1 kW (1.5 HP)
Maximum manometric head: 14.3m M.WC.



Laboratory Instrument India