



Email : sales@laboratoryinstrumentindia.com

Product Name :
Semi Auto Coagulation Analyzer

Product Code :
LBNY-0001-100002



Description :

Semi Auto Coagulation Analyzer

Technical Specification :

Ideal for low-volume testing or as a perfect back-up solution for fully automated optical systems requiring an alternative method for difficult samples, a reliable solution with an extensive test menu for all clotting assays from routine to specialized assays.

Parameters:

PT, APTT, Fibrinogen , Thrombin Time, Reptilase Time, Extrinsic pathway factors, Intrinsic pathway factors, Protein C; Protein S, Lupus Anticoagulants

Viscosity-Based (mechanical) Detection System: Immediate delivery of accurate and precise results for any type of coloured plasma, maximum precision for weak clot detection and standardization

Intuitive and standardised user interface within the Max Generation through a 7" color touch-screen

Management of Quality Control with Levey-Jennings graphs

Utility Requirement: Electrical power supply: 220VAC \pm 10%.

Suitable UPS with maintenance free batteries for minimum one-hour back-up

Laboratory Instrument India

