



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pascals Module

**Product Code :**  
LBNY-0005-10300024



**Description :**

**Structural Specifications:**

Quick connections for adaptation to feed hydraulics source.  
Stainless structure.  
Diagram in the front panel with similar distribution to the elements in the real unit.  
Screws, nuts, plates and all the metallic elements in stainless steel.

**Technical Specification :**

**Technical Data:**

Anodized aluminum profile frame and painted steel panel that guarantees good stability and resistance to the environment.  
Three vessels of 230 mm high:  
Inverted conical vessel with internal diameter from 26 mm. to 10 mm.  
Straight shaped vessel with internal diameter of 26 mm.  
Conical vessel with internal diameter from 26 mm. to 80 mm.



**Laboratory Instrument India**

