



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Depression Measurement System

**Product Code :**  
LBNY-0005-10300025



**Description :**

**Structural Specifications:**

Diagram in the front panel with similar distribution to the elements in the real unit.  
Quick connections for adaptation to feed hydraulics source.  
Stainless structure.  
Screws, nuts, plates and all the metallic elements in stainless steel.

**Technical Specification :**

**Technical Data:**

Anodized aluminum structure.  
Quick connections.  
Pressure-vacuum gauge adjusted from -1 to 0 bar.  
Dimensions: 220 x 110 x 420 mm. approx.  
Weight: 2 Kg.



**Laboratory Instrument India**

