



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Laminar Flow Demonstration

**Product Code :**  
LBNY-0005-10300026



**Description :**

Diagram in the front panel with similar distribution to the elements in the real unit.  
Quick connections for adaptation to feed hydraulics source.  
Stainless structure.  
Screws, nuts, plates and all the metallic elements in stainless steel.

**Technical Specification :**

**Technical Data:**

Depth of the table: adjustable depending on the models.  
Capacity of dye tank: 0.3 l.  
Width/length of the table: 400/210 mm.

**Hydrodynamic models:**

Two rectangular ones of 25 x 25 and 50 x 50 mm.  
Two circular ones of 25 and 50 mm. diameter.  
Wedge.  
Anodized aluminum structure.  
Easy and quick coupling system built-in.  
Dimensions: 870 x 450 x 400 mm. approx.  
Weight: 10 Kg.



**Laboratory Instrument India**