



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Water Hammer

**Product Code :**  
LBNY-0005-10300028



**Description :**

**Structural Specifications:**

Quick connections for adaptation to feed hydraulics source.  
Diagram in the front panel with similar distribution to the elements in the real unit.  
Stainless structure.  
Screws, nuts, plates and all the metallic elements in stainless steel.

**Technical Specification :**

**Technical Data:**

Connections system to the Hydraulics Bench or Basic Hydraulic Feed System with fast plugs.  
Pipe circuits in PVC.  
Valves to select the circuit.  
Constant level deposit, in methacrylate.  
Unload deposit, in methacrylate.  
2 adjustable equilibrium chimneys and subjection clips.  
Anodized aluminum structure.  
Easy and quick coupling system built-in.  
Dimensions: 1215 x 270 x 1430 mm.  
Weight: 15 Kg.

---

# Laboratory Instrument India

