



Email : sales@laboratoryinstrumentindia.com

Product Name :
Basic Hydraulic Feed System

Product Code :
LBNY-0005-10300030



Description :

Structural Specifications:

Screws, nuts, plates and all the metallic elements in stainless steel.
Diagram in the front panel with similar distribution to the elements in the real unit.
Stainless structure.
Quick connections for adaptation to feed hydraulics source.

Technical Specification :

Technical Data:

Stainless steel impeller.
Tank capacity: 140 l. approx.
Safety switch ON/OFF.
Supports for accommodating the test module.
This unit incorporates wheels for its mobility.
Flow meter.
Membrane type flow adjusting valve.
Centrifugal pump: 0.37 KW, 30-80 l./min at 20.1-12.8m., single-phase 220V./50Hz. or 110V./60Hz.
Dimensions: 1000 x 600 x 700 mm.
Weight: 40 Kg.

Laboratory Instrument India

