




Email : sales@laboratoryinstrumentindia.com

Product Name :
Flow Channel One M Length

Product Code :
LBNY-0005-10300032



Description :

Structural Specifications:

Diagram in the front panel with similar distribution to the elements in the real unit.
Quick connections for adaptation to feed hydraulics source.
Stainless structure.
Screws, nuts, plates and all the metallic elements in stainless steel.

Technical Specification :

Technical Data:

Channel of rectangular section with transparent walls in methacrylate, length: 1 m.
Tank with soothing of flow.
Easy and quick coupling system built-in.
Rigid and flexible pipes.
Anodized aluminum structure.
Regulating valves.
Wide range of available accessories.
Storage tank.
Dimensions: 1500 x 500 x 500 mm.
Weight: 40 Kg.

Laboratory Instrument India

