



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydraulic Power Supply

**Product Code :**  
LBNY-0005-10300039



www.laboratoryinstrumentindia.com

**Description :**

The Hydraulic Power Unit supplies hydraulic fluid to hydraulic learning systems including Basic Hydraulic System, Hydraulic Mechanic System, Electro-hydraulic Servo System, and the Combination Hydraulic Electric Servo System. This unit provides a minimum of 2.5 gpm flow at a maximum pressure to 500 PSI/ 3450 kPa.

**Features:**

Suction filter, strainer type  
Pressure Gauge, 2 in. size, 0-1000 PSI, liquid filled, with safety relief  
Electric motor starter, with start and stop push buttons  
Oil level gauge with temperature gauge  
Relief valve, pilot-operated, cartridge type  
Filler-Breather cap  
Power cord, grounded  
Return Hose  
Supply Hose  
Pump, fixed gear type  
Reservoir, 5 gal/19 Liters, welded square.

**Technical Specification :**

Hydraulic Power Supply



**Laboratory Instrument India**