



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Piping Learning System

**Product Code :**  
LBNY-0005-10300044



**Description :**

The Piping Learning System teaches how to layout and assemble fluid piping systems used in industrial, construction, agricultural, and transportation applications. The Piping Learning System Piping Technology learning system includes a mobile workstation, pump testing network, valve package, fitting package, fitting identification kits, Sloan valve maintenance kit, pipe insulation and hanger components, fitting storage, student learning materials for both theory and lab, and teachers guide. The workstation has slotted, movable members to enable a variety of hydraulic hose fluid piping application scenarios to be set up. Students will learn how to select, size, identify, and install a variety of types of industrial fluid piping, fittings and valves including: iron pipe, steel tubing, hydraulic hose, plastic pipe, copper tubing, globe valves, gate valves, check valves, and Sloan valves.

**Features:**

- Built-in Pump Leak Test System
- Built-in Storage
- Mobile Workstation
- Moveable Members.

**Technical Specification :**

Piping Learning System

---

# Laboratory Instrument India

