



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Compact Multi Pump Test Set Inverter Speed Control

**Product Code :**  
LBNY-0005-10300048



**Description :**

The unit consists of four different types of pump namely-centrifugal, turbine, reciprocating and axial flow pumps and a common storage tank, measuring instruments. Compact Multi Pump Test Set Inverter Speed Control is a compact self contained set for studying the characteristics for various types of pumps. By manipulating flow control valves and a selector switch, each pump can be operated individually. The unit is floor standing with removeable supports. Speed control is by an advanced inverter.

**Technical Specification :**

**Experiments:**

Pump input, output, and overall efficiency at various speeds.  
Flow rate vs head at various speeds.  
Effect of pulse chamber for reciprocating pump.



**Laboratory Instrument India**