



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Laboratory Washer and Dryer

**Product Code :**  
LBNY-0001-500017



**Description :**

Laboratory Washer and Dryer

**Technical Specification :**

Laboratory Washer And Dryer

General Features:

Application: Washing of chemical and microbiological glassware

Requirements

Under-bench, stainless steel, laboratory glassware washing machine, front loading with two tier system

Temperature on wash cycle and hot deionised water cycle 95oC

Temperature for drying cycle 110oC

Power supply : 220V-50 Hz or 415V 3 phase

Multi programmable (4), to include as a pre-wash, wash, water rinse, deionised water rinse (cold), deionised water rinse (hot) and drying cycles

Possibility to use adapted dispenser for liquid cleaning agents

Fitted with a steam condenser

To be supplied with racks and baskets allowing the ability to wash assorted sizes of flasks, beakers, cylinders and bottles.

Laboratory Washer and Dryer, Laboratory Washer and Dryer Equipments, Laboratory Washer and Dryer Tools, Laboratory Washer and Dryer Tool Kits, Laboratory Washer and Dryer Manufacturers, Laboratory Washer and Dryer Suppliers from India, China, Kenya



**Laboratory Instrument India**