



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifuge Tabletop

Product Code :
LBNY-0001-500018



Description :

Centrifuge Tabletop

Technical Specification :

Centrifuge Table Top

General Features:

Centrifuge for general chemistry use

Requirements

Capacity: 2 litres

Has digital display with actual speed and time

Variable speed and time control

Microprocessor controlled with automatic safety locks

Timer 0-60 min and setting for continuous operation.

Speed: 6000 rpm

Both timed and continuous operation

Safety device preventing the opening of the lid during centrifugation

In balance detector

To be supplied with

1 angle rotor for 8 tubes of 15/25 and 50 ml

Swing out rotor with 4 buckets

16 x 15ml; 8 x 25ml and 4 x 50 ml inserts

4 adapters/inserts for each 4 tubes of 15 ml for the swing out rotor bucket

Power supply: 220V-50 Hz.

Centrifuge Tabletop, Centrifuge Tabletop Equipments, Centrifuge Tabletop Tools, Centrifuge Tabletop Tool Kits, Centrifuge Tabletop Manufacturers, Centrifuge Tabletop Suppliers from India, China, Kenya



Laboratory Instrument India