



Email : sales@laboratoryinstrumentindia.com

Product Name :
Compact Centrifugal Pump Test Set Multi Speed

Product Code :
LBNY-0005-10300061



Description :

The unit consists of a storage tank and a centrifugal pump with multi speed motor on a steel base, and measuring instruments. Compact Centrifugal Pump Test Set Multi Speed Compact Centrifugal Pump Test Set Multi Speed is a self contained bench top unit for studying the centrifugal pump flow rate vs head at three different speeds.

Experiments:

Flow rate vs head at three different speeds.

Technical Specification :

Technical Data:

Pump ratings
Speeds : 3 different speeds.
Maximum head : 7 m water.
Maximum flow rate : 70 lpm.



Laboratory Instrument India

