



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydraulic Group

**Product Code :**  
LBNY-0005-10300068



www.laboratoryinstrumentindia.com

**Description :**

The aim of Hydraulic Group is to grant autonomy to other equipment included in Fluid Mechanics lab in order to have also the possibility to perform simultaneously different experiments without depending on the main hydraulic bench. Hydraulic Group designed to perform as an autonomous and portable unit for the supply of hydraulic power.

**Main Features:**

Storage tank with cap to prevent the inlet of impurities into the water.  
Easy operation.  
Pressure outlets for pump performance tests.  
Overflow chamber to avoid air inlet into the circuit.

**Technical Specification :**

**Technical Data:**

Tank:  
Capacity: 45 liters  
Material: Polypropylene  
Accessories: made of PVC

**Pump Features:**

Maximum manometric head: 23 T.M.H.  
Power consumption: 750W  
Max power: 950W

---

Motor speed: 2.900 rpm  
Volume 20 / 160 l/min  
H 21 / 10 T.M.H.  
H max: 23 T.M.H.  
H min: 10 T.M.H.



**Laboratory Instrument India**