




Email : sales@laboratoryinstrumentindia.com

Product Name :
Bernoullis Theorem Demonstration

Product Code :
LBNY-0005-10300077



Description :

Simple equipment suitable to study accurately the Bernoulli's Theorem Demonstration the system includes a manometer board with multiple tubes granting the possibility to read simultaneously different pressures along the pipe. The manometer has non leaking quick fittings in order to avoid the outlet of water when the system is disconnected from the hydraulic bench (not included in this item) is made through threaded screws and no tools are required.

Technical Specification :

Performable Experiments:

Demonstration of Bernoulli equation along the Venturi tube.
Study of static, dynamic and total pressures.
Pressure loss calculation in a Venturi tube.



Laboratory Instrument India