



Email : sales@laboratoryinstrumentindia.com

Product Name :
Roughened Bed

Product Code :
LBNY-0005-10300081



www.laboratoryinstrumentindia.com

Description :

The flow bed is made of stainless steel to which pyramid-shaped epoxy resin bodies are bonded. Using the Roughened Bed accessory, the user can simulate a very rough flume bed in the Modular Flow Channel. Removable carrying handles are used to transport the unit to the place of use; the bed is bolted to the base of the channel. These bodies absorb flow energy from the water. Inlet and outlet panels are fitted at the ends.

Technical Specification :

Bed made of stainless steel, with bonded epoxy resin bodies
Flow bed for inserting in the Modular Flow Channel
Lxwxh 2500x300x90mm
Removable transport aids
Dimensions and Weight
l x w x h: 2500 x 300 x 90 mm
Weight: 18 kg.



Laboratory Instrument India

