



Email : sales@laboratoryinstrumentindia.com

Product Name :
Plain Beach Simulation

Product Code :
LBNY-0005-10300084



www.laboratoryinstrumentindia.com

Description :

The stainless steel experimental arrangement forms a smooth, impermeable beach. The inclination of the beach body can be continuously adjusted. The experimental layout consists of a beach body, a mounting plate and an inclination adjustment bar with a locking bolt. The accessory is intended to enable the user to investigate the behaviour of waves on a rising beach; the unit is used in combination with the Wave Generator. In this way different steeply rising beaches are very easy to simulate.

Technical Specification :

Beach simulation for installation in the Modular Flow Channel
Inclination up to a maximum of 18 provided.
Unit made of stainless steel, l x w x h 1300x304x480mm.

Technical Data:

Angle of inclination: 18
Dimensions and Weight:
l x w x h: 1300 x 304 x 480 mm
Weight: 14 kg.

Laboratory Instrument India

