



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Set Of Piers Three Profiles

**Product Code :**  
LBNY-0005-10300089



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Using the Set of Piers 3 Profiles accessory, it is possible to perform experiments on pillars in rivers in the Modular Flow Channel. The experimental set-up is made of PVC and is attached to the base of the channel using bolts. Pointed, rounded and flat flow splitters are provided; these can be combined as required. The unit contains 3 different, separable pillar types, 2 of each type of pillar fitted to the base plate represent the complete experimental set-up ready for installation in the flow channel.

**Technical Specification :**

Flow splitters for installing in the Modular Flow Channel  
Sets of flow splitters: pointed, rounded, flat, l x w x h 100x30x300mm  
Experimental set-up made of PVC, l x w x h 330x160x315mm.

**Technical Data:**

Rounded flow splitter  
Radius: 15mm  
Pointed flow splitter  
Included angle: 41  
Dimensions and Weight  
l x w x h: 330 x 160 x 315 mm  
Weight: 5 kg.



**Laboratory Instrument India**