



Email : sales@laboratoryinstrumentindia.com

Product Name :
Set Of Piers Three Profiles

Product Code :
LBNY-0005-10300089



www.laboratoryinstrumentindia.com

Description :

Using the Set of Piers 3 Profiles accessory, it is possible to perform experiments on pillars in rivers in the Modular Flow Channel. The experimental set-up is made of PVC and is attached to the base of the channel using bolts. Pointed, rounded and flat flow splitters are provided; these can be combined as required. The unit contains 3 different, separable pillar types, 2 of each type of pillar fitted to the base plate represent the complete experimental set-up ready for installation in the flow channel.

Technical Specification :

Flow splitters for installing in the Modular Flow Channel
Sets of flow splitters: pointed, rounded, flat, l x w x h 100x30x300mm
Experimental set-up made of PVC, l x w x h 330x160x315mm.

Technical Data:

Rounded flow splitter
Radius: 15mm
Pointed flow splitter
Included angle: 41
Dimensions and Weight
l x w x h: 330 x 160 x 315 mm
Weight: 5 kg.



Laboratory Instrument India