



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Fittings

Product Code :
LBNY-0005-10300097



Description :

The models supplied are made of plastic and are equipped with connections for pressure measurement; these are connected to the Computer Linked Air Flow Bench. It is possible to investigate the pipe friction losses in a flow of air for different fittings such as a straight pipe section, an elbow, a bend, and different pipe inlets.

Technical Specification :

Straight pipe section, extension possible
Two different pipe inlets
Model on pipe friction losses
Lxwxh unassembled 1500x200x400mm,5kg
Bend and elbow, each 90.

Technical Data:

Pipe inlet: undisturbed inlet 22mm radius
Discontinuous enlargement: 35...53.6mm
Pipe section with flange:
1235mm long, inside diameter 53.6mm
Extension: 991mm long (with inlet)
Dimensions and Weight of Pipe Fittings
Lxwxh unassembled: 1500 x 200 x 400 mm
Weight: 5 kg.



Laboratory Instrument India