



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pipe Fittings

**Product Code :**  
LBNY-0005-10300097



**Description :**

The models supplied are made of plastic and are equipped with connections for pressure measurement; these are connected to the Computer Linked Air Flow Bench. It is possible to investigate the pipe friction losses in a flow of air for different fittings such as a straight pipe section, an elbow, a bend, and different pipe inlets.

**Technical Specification :**

Straight pipe section, extension possible  
Two different pipe inlets  
Model on pipe friction losses  
Lxwxh unassembled 1500x200x400mm,5kg  
Bend and elbow, each 90.

**Technical Data:**

Pipe inlet: undisturbed inlet 22mm radius  
Discontinuous enlargement: 35...53.6mm  
Pipe section with flange:  
1235mm long, inside diameter 53.6mm  
Extension: 991mm long (with inlet)  
Dimensions and Weight of Pipe Fittings  
Lxwxh unassembled: 1500 x 200 x 400 mm  
Weight: 5 kg.



**Laboratory Instrument India**