



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flow Around A Bend

Product Code :
LBNY-0005-10300098



Description :

A transparent pipe bend with a constant square cross-section is attached to the test stand using quick action fasteners. The measuring glands are connected to the Multi Tube Manometer. The pipe bend is equipped with measuring glands at regular intervals for measuring the pressure loss.

Technical Specification :

Pipe bend made of Plexiglass, with pressure measuring glands
Experimental set-up for determining the pressure loss in a 90 pipe bend through which there is a flow
Accessory for the Air Flow Bench.

Technical Data:

Number of measuring glands: 28
Pipe bend cross-section: 50x100mm
Dimensions and Weight of Flow Around a Bend
l x w x h : 230 x 180 x 200 mm
Weight: 4 kg.



Laboratory Instrument India

