



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pump Test Accessory

Product Code :
LBNY-0005-103000102



Description :

The unit incorporates a pressure sensor to measure the inlet condition of the water. The consists of a return tube incorporating a gate valve, which can be attached to the pump outlet on the service unit.

Technical Specification :

Mounts on a dedicated service unit
Compact accessory to enable investigation of pump performance
Links to a PC via the service unit and a USB interface console

Experimental Capabilities:

Determination of characteristic performance curves for a peripheral pump, including constant speed head/flow and efficiency curves.



Laboratory Instrument India