



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Reciprocating Compressor Module

**Product Code :**  
LBNY-0005-103000103



**Description :**

**Operating Conditions:**

Operating relative humidity range: 80% at temperatures Operating environment: Laboratory  
Operating temperature range: +5C to +40C  
Storage temperature range: 25C to +55C (when packed for transport).

**Technical Specification :**

**Experiments:**

Energy balance for a compressor  
Mechanical, volumetric and isothermal efficiencies  
Variation of compressor performance with pressure  
Variation of compressor performance with speed  
Thermodynamics of a compressor.



**Laboratory Instrument India**