



Email : sales@laboratoryinstrumentindia.com

Product Name :
Reciprocating Compressor Module

Product Code :
LBNY-0005-103000103



Description :

Operating Conditions:

Operating relative humidity range: 80% at temperatures Operating environment: Laboratory
Operating temperature range: +5C to +40C
Storage temperature range: 25C to +55C (when packed for transport).

Technical Specification :

Experiments:

Energy balance for a compressor
Mechanical, volumetric and isothermal efficiencies
Variation of compressor performance with pressure
Variation of compressor performance with speed
Thermodynamics of a compressor.



Laboratory Instrument India