



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Positive Displacement Pump Support Module

**Product Code :**  
LBNY-0005-103000105



**Description :**

**Operating Conditions:**

Operating environment: Laboratory

Operating relative humidity range: 80% at temperatures Operating temperature range: +5C to +40C

Storage temperature range: 25C to +55C (when packed for transport).

**Technical Specification :**

**Technical Specifications:**

Oil type: Shell Morlina 10 or equivalent

Oil reservoir capacity: 30 litres

Temperature measurement: K-type thermocouple and display

Nett dimensions and weight (assembled): 1540 mm x 1620 mm x 660 mm and 138 kg (with no oil in the reservoir)

Approximate packed volume and weight: 1.6 m3 and 250 kg (with oil).



**Laboratory Instrument India**

