



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Centrifugal Fan Module

**Product Code :**  
LBNY-0005-103000106



**Description :**

**Experiments:**

Determination of the specific speed of the fan  
Performance of a centrifugal fan  
Variation of fan performance with type of impeller  
Variation of fan performance with speed  
Non-dimensional performance curves.

**Technical Specification :**

**Technical Specifications:**

Operating environment: Laboratory  
Operating relative humidity range: 80% at temperatures  
Operating temperature range: +5C to +40C  
Storage temperature range: 25C to +55C (when packed for transport).



**Laboratory Instrument India**