



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cavitation Demonstration Unit

Product Code :
LBNY-0005-103000107



Description :

Operating Conditions:

Operating environment:Laboratory

Operating relative humidity range: 80% at temperatures Operating temperature range:+5C to +40C

Storage temperature range:25C to +55C (when packed for transport).

Technical Specification :

Experiments:

Demonstrations of cavitation

Investigations into cavitation and the Venturi, including: Flow and pressure in the Venturi

Essential Services

How to predict the onset of cavitation

Floor space needed:Approximately 1 m x 1.5 m of solid, level floor

Electrical supply: Single-phase earthed electrical supply, 230 VAC, 50 Hz.



Laboratory Instrument India

