



Email : sales@laboratoryinstrumentindia.com

Product Name :
Reynolds No And Transitional Flow

Product Code :
LBNY-0005-103000111



Description :

Experiments:

Investigation of the effect of varying viscosity and demonstration that the Reynolds number at transition is independent of viscosity.
Determination of transition Reynolds numbers and comparison with accepted values.
Demonstration of transition between laminar and turbulent flow.

Technical Specification :

Operating Conditions:

Operating environment: Laboratory
Operating relative humidity range: 80% at temperatures
Operating temperature range: +5C to +40C
Storage temperature range: 25C to +55C (when packed for transport).



Laboratory Instrument India