



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Series And Parallel Pump Test Set

**Product Code :**  
LBNY-0005-103000114



**Description :**

**Experiments:**

Parallel operating characteristics of two pumps operating at different speeds  
Series operating characteristics of two pumps operating at different speeds  
Performance of single centrifugal pump  
Series operating characteristics of two similar pumps  
Parallel operating characteristics of two similar pumps.

**Essential Base Unit:**

Volumetric Hydraulic Bench or  
Gravi metric Hydraulic Bench.

**Technical Specification :**

**Essential Services:**

Water supply: From the Hydraulic Bench  
Floor space needed: 1 m x 1 m (plus space for the hydraulic bench)  
Electrical supply: 230 VAC, 50/60 Hz, 2 A, single-phase or two-phase  
Operating Conditions  
Operating environment: Laboratory  
Storage temperature range: 25C to +55C (when packed for transport).



**Laboratory Instrument India**