



Email : sales@laboratoryinstrumentindia.com

Product Name :
Series And Parallel Pump Test Set

Product Code :
LBNY-0005-103000114



Description :

Experiments:

Parallel operating characteristics of two pumps operating at different speeds
Series operating characteristics of two pumps operating at different speeds
Performance of single centrifugal pump
Series operating characteristics of two similar pumps
Parallel operating characteristics of two similar pumps.

Essential Base Unit:

Volumetric Hydraulic Bench or
Gravi metric Hydraulic Bench.

Technical Specification :

Essential Services:

Water supply: From the Hydraulic Bench
Floor space needed: 1 m x 1 m (plus space for the hydraulic bench)
Electrical supply: 230 VAC, 50/60 Hz, 2 A, single-phase or two-phase
Operating Conditions
Operating environment: Laboratory
Storage temperature range: 25C to +55C (when packed for transport).



Laboratory Instrument India