



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Axial Flow Pump Apparatus

**Product Code :**  
LBNY-0005-103000116



**Description :**

Determination of the flow rate with the electromagnetic flow rate sensor  
Water tank with sight glass  
Digital display of flow rate, pressure and temperature  
Investigation of an axial-flow pump  
Powered by Universal Drive and Brake Unit  
Closed water circuit.

**Technical Specification :**

**Technical Data:**

Pump  
Max. head: 1,75m  
Max. flow rate: 700L/min  
Output: 1000W at 1000min-1  
Tank: 160L  
Measuring ranges  
Pressure (inlet): -1...1bar  
Pressure (outlet): 0...0,6bar  
Flow rate: 0...833L/min  
Temperature: 0...100C.



**Laboratory Instrument India**