



Email : sales@laboratoryinstrumentindia.com

Product Name :
Axial Flow Pump Apparatus

Product Code :
LBNY-0005-103000116



Description :

Determination of the flow rate with the electromagnetic flow rate sensor
Water tank with sight glass
Digital display of flow rate, pressure and temperature
Investigation of an axial-flow pump
Powered by Universal Drive and Brake Unit
Closed water circuit.

Technical Specification :

Technical Data:

Pump
Max. head: 1,75m
Max. flow rate: 700L/min
Output: 1000W at 1000min-1
Tank: 160L
Measuring ranges
Pressure (inlet): -1...1bar
Pressure (outlet): 0...0,6bar
Flow rate: 0...833L/min
Temperature: 0...100C.



Laboratory Instrument India