



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cooling Tower 115kW

Product Code :
LBNY-0005-103000117



Description :

Water flow measurement at the inlet and outlet
Combined temperature/humidity measurement
Wet cooling tower with fan and pump for operation
Use of corrosion-resistant materials.

Technical Specification :

Technical Data:

Cooling tower
Water flow rate 9m/h
Cooling capacity approx. 115kW
Pump max.: 19,2m/h
Fan max.: 3900m/h at 1450min-1
Measuring ranges;
Humidity: 2x 5...95% r.h.
Pressure: 1x 0...2,5bar, 1x 0...6bar
Temperature: 2x -10...50C, 2x 0...80C
Flow rate: 0...20m/h.

Laboratory Instrument India

