



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vapour Pressure Of Water Marcet Boiler

**Product Code :**  
LBNY-0005-103000125



**Description :**

Measuring a vapour pressure curve for saturated vapour  
Bourdon tube pressure gauge to indicate pressure  
Digital temperature display  
Boiler with insulating jacket  
Temperature limiter and safety valve protect against overpressure in the system.

**Technical Specification :**

**Technical Data:**  
Bourdon tube pressure gauge: -1...24bar  
Heater: 2kW  
Safety valve: 20bar  
Temperature limiter: 200C  
Boiler, stainless steel: 2L  
Measuring ranges  
Pressure: 0...20bar  
Temperature: 0...250C.

---

# Laboratory Instrument India

